


# Guide to needle size

Use the shortest thinnest needle you can



LENGTH	<b>30G YELLOW</b> ∅   0.3 mm diameter	<b>29G RED</b> ∅   0.33 mm diameter	<b>27G GREY</b> ∅   0.4 mm diameter	<b>26G BROWN</b> ∅   0.45 mm diameter	<b>25G ORANGE</b> ∅   0.5 mm diameter	<b>24G PURPLE</b> ∅   0.55 mm diameter	<b>23G BLUE</b> ∅   0.6 mm diameter	<b>22G BLACK</b> ∅   0.7 mm diameter	<b>21G GREEN</b> ∅   0.8 mm diameter	<b>18G PINK</b> * DRAWING UP ONLY *	LENGTH
12 mm	12mm [½ inch] Nevershare/Unisharp	12mm [½ inch] Unisharp Fixed	12mm [½ inch]	12mm [½ inch]							12 mm
16 mm			16mm [⅝ inch]		16mm [⅝ inch] Total Dose		16mm [⅝ inch]				16 mm
20 mm			20mm [¾ inch]								20 mm
25 mm				25mm [1 inch]	25mm [1 inch] Total Dose	25mm [1 inch] Unisharp	25mm [1 inch] Unisharp				25 mm
30 mm							30mm [1¼ inch] Total Dose	30mm [1¼ inch]			30 mm
40 mm			40mm [1½ inch]		40mm [1½ inch]			40mm [1½ inch]	40mm [1½ inch] Total Dose	40mm [1½ inch] blunt/sharp	40 mm
50 mm									50mm [2 inch]		50 mm