


Guide to needle size

Use the shortest thinnest needle you can



LENGTH	30G YELLOW ∅ 0.3 mm diameter	29G RED ∅ 0.33 mm diameter	27G GREY ∅ 0.4 mm diameter	26G BROWN ∅ 0.45 mm diameter	25G ORANGE ∅ 0.5 mm diameter	24G PURPLE ∅ 0.55 mm diameter	23G BLUE ∅ 0.6 mm diameter	22G BLACK ∅ 0.7 mm diameter	21G GREEN ∅ 0.8 mm diameter	18G PINK * DRAWING UP ONLY *	LENGTH
12 mm	12mm [1/2 inch] Unisharp & Unisharp Fixed	12mm [1/2 inch] Unisharp Fixed	12mm [1/2 inch] Unisharp & Unisharp Fixed	12mm [1/2 inch] Unisharp							12 mm
16 mm			16mm [5/8 inch]		16mm [5/8 inch] Total Dose		16mm [5/8 inch]				16 mm
20 mm			20mm [3/4 inch]								20 mm
25 mm				25mm [1 inch] Unisharp	25mm [1 inch] Total Dose	25mm [1 inch] Unisharp	25mm [1 inch] Total Dose		25mm [1 inch] Unisharp		25 mm
30 mm							30mm [1 1/4 inch] Total Dose	30mm [1 1/4 inch]			30 mm
40 mm			40mm [1 1/2 inch]		40mm [1 1/2 inch] Unisharp			40mm [1 1/2 inch]	40mm [1 1/2 inch] Total Dose	40mm [1 1/2 inch] blunt/sharp	40 mm
50 mm									50mm [2 inch]		50 mm